

# Stop snoring. Sleep better. Get aveoTSD®.

**aveoTSD®**  
Tongue Stabilizing Device

**The brilliantly  
simple solution for  
problem snoring.**




**aveoTSD®**

**How aveoTSD works**



TONGUE




In this MRI\* image, the tongue falls into the back of the airway as a person sleeps. This blocks the airway, leading to snoring.

\*Magnetic Resonance Imaging  
GE Signa Profile EXCITE 0.2T

BLOCKED  
AIRWAY

**aveoTSD®**  
Tongue Stabilizing Device



This MRI image shows the aveoTSD holding the tongue gently forward. Note how the airway is now open and clear. This stops or greatly reduces snoring.

OPEN  
AIRWAY

### **Patients:**

*Please print this sheet and take it to your dentist or physician who can prescribe you an aveoTSD.*

### **Dentists/Physicians:**

*Call 800-334-1979 to order.  
Visit [www.getaveo.com](http://www.getaveo.com) for more info.*

**800-334-1979 • [www.getaveo.com](http://www.getaveo.com)**