

## Revolutionary PFM Strength & Esthetics




Obsidian® Pressed to Metal – non-precious: \$139/unit\*

- FIVE TIMES STRONGER than Ceramco®3 PFM ceramic<sup>†</sup>
- TWICE THE CHIP RESISTANCE of Ceramco3 PFM ceramic<sup>†</sup>
- **IDEAL FOR BRIDGES**, cases with insufficient vertical space for full-contour zirconia, and screw-retained restorations











Chip resistance is measured by applying a point-force (load) on the edge of a test sample until chipping occurs. Obsidian Pressed to Metal withstands more than twice the load compared to feldspathic porcelains.





## "The best looking PFMs I've seen in 27 years as a ceramist."

- Joe Hattouni, CDT, General Manager Fixed Prosthodontics







After



Before



After



Obsidian
30-month Clinical Performance
+ + + + +



## LONG-TERM CLINICAL PERFORMANCE

The performance and esthetics of *Obsidian* crowns during the 30-month evaluation period has been exceptional. The restorations will continue to be monitored over time. *Obsidian* received a clinical performance rating of 98% at the 30-month recall.

## Sledgehammer Toughness Test at www.obsidianceramic.com



Obsidian® (left), Ceramco®3 (right)



Crowns hit with sledgehammer



Obsidian (left), Ceramco3 (right)

\*Price does not include shipping or applicable taxes. †Data provided by Glidewell Laboratories Research and Development Department, 2015.

Obsidian is a registered trademark of Prismatik Dentalcraft, Inc. Ceramco is a registered trademark of DENTSPLY Ceramco.

For more information
800-854-7256
www.glidewelldental.com

